



Patient Status Engine

AI based on Vital Signs

*StartUp- und Digitalisierungspreis der
ENTSCHEIDERFABRIK*

Medica 2018

Digital Health?

Digitisation of Process and Information

- Hospital information systems (Electronic patient records etc)
- Electronic prescribing, billing, imaging etc
- Many apps – reminders for appointments, medications etc

Combine with connectivity

Business process engineering




Digital Health?

“From ineffective electronic health records, to an explosion of direct-to-consumer digital health products, to apps of mixed quality.

This is the digital snake oil of the early 21st century.”

Dr James Madara

American Medical Association CEO



This clock has stopped.




But it is extraordinarily accurate. How?

What do these numbers represent?


98.6

37.0



These numbers?

120 / 80



According to the blood pressure classification by the WHO/ISH.

These numbers?

60 - 80



Are these numbers Right or Wrong?


Sometimes these numbers are right

Most of the time they are wrong
and often misleading


It Depends on




Time of Day (Circadian Cycles)




Age



Activity / Exercise (and where it's measured – *Ouch!*)




Time of Month (Menstrual Cycles)




© Microsoft Corporation. All Rights Reserved.

And Not Only Temperature




Our bodies have more than 100 circadian rhythms

Looking Continuously




And Looking More Closely




Reveal What's Really Going On

A: Mean 65.0 SD = 4.8



Healthy

B: Mean 65.0 SD = 4.7



Sleep Apnoea

Digital Health?

Digitisation of Content

- Genome
- Epigenome
- What about your Physionome?

The truly disruptive step






The Digital Future of Patient Monitoring

- *Anytime, anywhere, in hospital, at home*
- *Robust, accurate, wireless patient data capture and analysis*
- *Direct, real-time, quantitative patient outcome measures*

Enabling the Solution: The Patient Status Engine



New Possibilities In Patient Monitoring

The Isansys Patient Status Engine PSE



Status quo:

- selective values
- time pressure
- mistakes & errors
- wired patients
- demobilizing
- interruptions






continuous

algorithm

timesaving

unburden

precise



location-independent

mobilizing

Transmission via
WLAN, GSM



The Patient Status Engine (PSE)


New Possibilities in Patient Monitoring

A complete system:



- **Wireless wearable sensors:**
 - *Heartrate, Temperature, Respiration Rate*
 - *SPO2 – Oxygen Saturation*
 - *Blood Pressure (as required)*
- **Documentation:**
 - *Automatically filled patient charts*
 - *fast, accurately and correctly*
- **Interoperable:**
 - *Integration into the digital hospital IT infrastructure / electronic patient records*
 - *Telemedicine: authorized access from everywhere, remotely*
- **Early Warning Scores (EWS)**
 - *Critical data at a glance*



Notification Score: Digitizing the Patient: Digital Clinical Decision Support



Notification Score: Single Vital Parameter




Current situation: Step 1 in digitizing the patient

- Advantage towards snapshots all 4 hours
- Physician has the chance to see patient status change
- Medical staff has still to be active with their “Experience Algorithm”



Notification Score: Single Parameter enhanced by Context




A further step:

- A combination of several vital data and additional information enables PSE to make more reliable estimation of actual patient condition



Notification Score: Single Vital Parameter over Course of Time




Next step towards digitization:

- Personal "Experience Algorithm" will be implemented as ward or hospital wide Standard Operating Procedure (SOP)





Notification Score: Single Vital Parameter over Course of Time

Our way to using Artificial Intelligence



Notification Score: Neuronal Network (Algorithm) for Patient Status





Next step towards digitization:

- Mapping complex input of vital parameters to establish standardized scores

Robust outcome in presence of noise and lack of input parameter required

Notification Score: Patient Status Engine's next evolutionary step

Neuronal Network (Algorithm) for Status Prediction



Digitizing the Patient: Digital Clinical Decision Support

- Medical experience of physicians is used through the digital health solution
- Algorithm can be used as SOP across clinical institutions
- New methods of medical learning possible



Contact

Isansys Lifecare Europe GmbH
Michael Heinlein
Charlottenstraße 72
14467 Potsdam
michael.heinlein@isansys.eu

Info.eu@isansys.com
<http://www.isansys.com>